

- Media Bottle Available in 5 different sizes: 60mL, 125mL, 250mL, 500mL, 1000mL
- Made entirely of 100% USP Class VI Components
- BPA Free PETG
- Compactable with 38-430 VersaCap® that are leakproof and provides flexibility of interchangeable connectors.
- Can withstand low temperatures of up to -70°C
- Manufactured in a Class 7 cleanroom.
- Dnase, Rnase & Pyrogen Free
- Assembled and Packaged in our ISO 13485 Facility in India

Foxx Life Sciences New **EZBio® Titanium Square Bottle** has enhanced the media bottle design and offers better lead time and pricing. The New Titanium Square Media Bottles are perfect for storing and sampling active pharmaceutical ingredients and bulk intermediates, preparing and storing buffers, and for culture. These bottles can also be used for solutions or long-term storage of pH-sensitive liquids. The **EZBio® Titanium Square Bottle** space saving configuration allows for maximum use of precious incubator space.

The **EZBio® Titanium Square Bottle** is molded from virgin, optically clear Polyethylene Terephthalate Glycol (PETG) BPA Free Resin. The **EZBio® Titanium Square Sterile Bottles** are certified as non-pyrogenic and DNase/RNase-free, ensuring they are free of contaminants. Easy to read raised graduation marks help in dispensing media.

Available storage sizes are 125/250/500/1000mL. Each bottle has a 38-430 neck compatible with our patented 38-430 **VersaCap®** assemblies and other 38-430 Caps standards. The **VersaCap®** Technology provides flexibility for choosing a wide range of port combinations use for mixing sampling or fluid transfer to closed systems. Foxx also offers the ability to provide and produce custom components, assemblies, and variations on time. For more information, please contact us at sales@foxxlifesciences.com or 1-603-890-3699.

- 1 38-430 VersaCap® Technology
- 2 Ultra-pure USP Class VI resins
- 3 Easy-to-read metric graduation marks certified to ±5% accuracy
- 4 Available storage sizes are 125, 250, 500, 1000mL



 BPA-Free

Non-pyrogenic 

USP
Class VI

 GAMMA
STERILIZED

DNasefree
RNasefree



ISO 13485



Part Number	Description	Material	Cap	Bottle Size	Cap Size (mm)	Sterile	Qty/Pack
627-1101-FLS	EZBio® Titanium Square Media Bottles	PETG	Yes	60mL	24-415	Yes	192/CS
627-2101-FLS	EZBio® Titanium Square Media Bottles	PETG	Yes	125mL	38-430	Yes	96/CS
627-3101-FLS	EZBio® Titanium Square Media Bottles	PETG	Yes	250mL	38-430	Yes	60/CS
627-4101-FLS	EZBio® Titanium Square Media Bottles	PETG	Yes	500mL	38-430	Yes	48/CS
627-5101-FLS	EZBio® Titanium Square Media Bottles	PETG	Yes	1000mL	38-430	Yes	24/CS
627-1003-FLS	EZBio® Titanium Square Media Bottles	PETG	No	60mL	24-415	No	192/CS
627-2003-FLS	EZBio® Titanium Square Media Bottles	PETG	No	125mL	38-430	No	96/CS
627-3003-FLS	EZBio® Titanium Square Media Bottles	PETG	No	250mL	38-430	No	60/CS
627-4003-FLS	EZBio® Titanium Square Media Bottles	PETG	No	500mL	38-430	No	40/CS
627-5003-FLS	EZBio® Titanium Square Media Bottles	PETG	No	1000mL	38-430	No	24/CS

