

- **EZLabpure™** Aspirator Bottles are manufactured using high purity virgin USP Class VI Medical grade polypropylene

- Very strong and durable for extended user over time. They are equipped with inner-cap for a tight seal. Ideal for containing liquid powders and granulates. Comes with high quality & durable stopcock or tap. These are used for Packing, Shipping & Storage purpose

- All can be used at a temperature below 0°C

- Aspirator Bottles are autoclavable at 121°C for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic.

Caution

- Do not store strong oxidizing agents
- Do not place any plastic labware in a direct flame
- Do not mix any chemicals that may result in a thermic reaction, which can cause product failure
- Consult local fire codes prior to storage of flammable liquids in carboys/bottles/aspirators
- Consult OSHA guidelines prior to handling or lifting filled carboys/bottles/aspirators

Foxx Life Sciences offers a premium line of high quality, polypropylene aspirator bottles. **EZLabpure™** Aspirator Bottles come with strong easy-grip handle and Spigot. Our versatile bottles are leak-proof and aim to take full advantage of storage efficiency at minimal effort. Overall functionality is available for many applications, in five liter, ten liter, and twenty liter sizes, with spigots and PP cap

Autoclavability	Yes
Translucent	Yes
Resistant	To be effected by radiations
Sterilization	Gamma
Leak Proof	Yes

- 1 PP (polypropylene) cap
- 2 Ergonomic handle design to make handling, lifting, and pouring bottle easier
- 3 Important instructions and warnings for the **EZLabpure™** Aspirator Bottle
- 4 High quality and durable stopcock with ON and OFF engraving



MADE IN  ISO 13485 BPA-Free 

USP
Class VI

FDA

AUTOCLAVABLE



EZLabpure™ Aspirator Bottles Ordering Information:

Product Number	Description	Capacity (L)	Mouth Diameter (mm)	Height (mm)	Autoclavable
155-J897-FLS	5L PP Aspirator Bottle with White Cap and Spigot	5	51.4	374	Yes
155-K897-FLS	10L PP Aspirator Bottle with White Cap and Spigot	10	51.4	445	Yes
155-L897-FLS	20L PP Aspirator Bottle with White Cap and Spigot	20	51.4	525	Yes

Disclaimer: Autoclave Carboy and Cap Separately

