



- **EZLabpure™** Round Carboys are manufactured using high purity virgin USP Class VI Medical grade polypropylene
- These carboys come with secure sealing cap for better security and long-term storage
- These are very strong and durable for extended use over time with a strong carrying handle on both sides of the carboy
- Comes with wide mouth easy withdrawal and refill of liquids, chemicals, reagents, culture media and standards
- PP carboys are autoclavable at 121°C, 15psi pressure, 20 mins.
- Can be used at a temperature below 0°C
- Superior chemical resistance

Caution

- Autoclave carboy and cap separately

MADE IN  **ISO 13485 BPA-Free**



AUTOCCLAVABLE

Foxx Life Sciences offers a premium line of high quality, round wide mouth polypropylene (PP) laboratory carboys. Rugged and durable, **EZLabpure™** Round Carboys are the safe, shatterproof alternative to glass for mixing, storing and dispensing liquids in volumes ten and twenty liters.

Product Number	Description
155-B187-FLS	10L Wide Mouth PP Carboy with White Cap
155-C187-FLS	20L Wide Mouth PP Carboy with White Cap

- 1 Wide 100mm mouth for use with solids, powders and overhead mixing equipment
- 2 Handle design to make handling, lifting and pouring carboy easier
- 3 Round carboy design maximizes mixing efficiency
- 4 Graduation marks available on 10L and 20L size





EZLabpure™ Round Carboys Also Available:

Foxx Life Sciences offers a premium line of high quality, round polypropylene (PP) Laboratory carboys with the patented VersaCap® technology. They can be used for mixing, storing and dispensing liquids in volumes from ten to twenty liters. These carboys are available with and without a spigot.

EZLabpure™ Round Carboys are compatible with Foxx's patented VersaCap® that is leak-proof, has replaceable platinum-cured silicone gaskets and provides the ultimate flexibility. The open VersaCap® allows a variety of configurations by changing the adapter insert, enabling users to quickly change the top connections of each carboy or bottle to match their application needs.



EZLabpure™ Wide Mouth Carboy Ordering Information:

Product Number	Description	Capacity (L)	Cap Size (mm)	Autoclavable	Qty (EA)
155-B117-FLS	10L Round PP Carboy with 83B Cap	10	83	Yes	1
155-C117-FLS	20L Round PP Carboy with 83B Cap	20	83	Yes	1
155-B817-FLS	10L Round PP Carboy with 83B Cap and Spigot	10	83	Yes	1
155-C817-FLS	20L Round PP Carboy with 83B Cap and Spigot	20	83	Yes	1